



Specification

Product Name : 1-Butyl-2,3-dimethylimidazolium
trifluoromethanesulfonate, 98%

Catalog No. : 111199

CAS No. : 765910-73-4

Molecular Formula : C₁₀H₁₇F₃N₂O₃S

Molecular Weight : 302.31

| Test | Standard |
|-------------|------------------------------|
| Appearance | : White to off-white solid |
| IR | : Conforms to structure |
| Purity | : ≥98% (with total nitrogen) |
| HPLC | : ≥98% |